



for  
Freeze Dryer Controlling



[ OPERATION ]  
[ MANUAL ]





# for Freeze Dryer Controlling

## OVERVIEW

The purpose of this document is to describe the purpose, implementation and use of the **FYRA Freeze Dryer Controller**.


The **FYRA for Freeze Drying Bundle** includes all of the components necessary to enable easy retrofits of older or broken freeze dryers. It also includes a vacuum control facility that can enhance drying time and improve product quality in many instances. This vacuum control ability is found typically in the most expensive freeze drying products.

## INCLUDES

- FYRA platform
- Temperature card to read type “K” thermocouples for displaying 2 temperatures: *Shelf* and *Ambient*
- Vacuum gauge driver card and Agilent 536 vacuum gauge
- VLC card and VacStable valve for controlling the vacuum pressure
- Solid state relay driver card
- UI that enables manual freeze-dry or automated freeze-drying modes
- **Variables** that can be customized to your drying process
- Platform designed from the ground up to be safe and reliable.
- Solid-state quad-relay

## INSTALLATION CONSIDERATIONS

- Ideally the bleed valve should be installed on a vessel fitting **opposite** where the vacuum pump is attached. This isn't required, but helpful to promote flow.
- If the freeze-dryer does not have Type K Thermocouples, Type K Thermocouples should be installed on the shelf.
  - For the shelf temperature, it is important that the thermocouple is *on the actual shelf* otherwise inconsistent heating and poor temperature control will result.
  - Type K Thermocouples should be placed where ambient temperature can be taken.



## VARIABLES

### For Manual and Auto Mode


- **Shelf Limit:** When in step 3, or *Drying mode*, this is the temperature that the shelf will be controlled at 70 -170 F
  - User settable
  - Common range for this is 20-40F, but can be just about anything
- **Vac Setpoint:** This is the vacuum level, expressed in Torr that the system will not go below.
  - Typical values for this range from 0.3 Torr to 1 Torr.
  - Vacuum level is maintained by bleeding in a small amount of air.
  - It is best that this valve be placed on the chamber opposite where the vacuum pump is installed.

### For Auto Mode Only

- **Freeze Time:** The time the user would like the freeze-dryer to spend freezing the product.
  - Note: The vacuum pump will be turned on the *last 60 minutes* of the user specified freeze time.
- **Dry Time:** The time the user would like the freeze-dryer to *Dry* the product. During *Drying*, the refrigeration system, the vacuum pump and the shelf heater are on.

## SAFETY

If temperature of shelf exceeds preset *Shelf Limit*, which might happen if a thermocouple fails or is disconnected, the heater will be turned off. The pump and refrigeration would remain on.



## 4 MODES

- **Freezing**
  - Refrigeration is active
- **Freezing + Pumping**
  - Both refrigeration system and vacuum pump are energized
- **Drying**
  - Refrigeration, vacuum pump and heater are on, so long as the heat doesn't exceed the shelf temp set point.
- **Defrost**
  - Heater is on, Temperature is factory set to 220 F



## AUTO

1. Turns refrigeration system on for "Freeze Time"
2. In last 60 minutes of Freeze Time, the vacuum pump turns on.
  - a. If Freeze Time < 60 minutes, then the vacuum pump and the refrigeration system will come on immediately.
3. After Freeze Time is complete, **Final Dry** will last for the time set in "Final Dry."
4. When Final Dry is on, the refrigeration system, the pump and the heater are on
5. The heater will stay at "shelf temp"
  - a. Heat *does not* turn on unless the shelf temperature is *below* the shelf temperature limit variable
6. Refrigeration pump and heat is off after the auto cycle is complete

## MANUAL

1. Touch screen has 4 buttons describing above modes [1. **Freezing**, 2. **Freezing + Pumping**, 3. **Drying**, and **Defrost**]
2. Modes can be manually switched between/changed by tapping on these buttons

